

ScanAfia XS and XTK Hospital Bed

NEW! The safe choice for hospital wards

The quality ScanAfia XS and XTK hospital bed is the safe choice for hospital wards. The bed has been designed for use by healthcare professionals in hospitals, healthcare centres and similar facilities, and it fulfils the requirements of the new EN 60601-2-52:2010 standard for medical beds.

Although it can take loads up to 260 kg, the ScanAfia XS and XTK hospital bed is still light and easy to manoeuvre, thanks to its double castors. The design places emphasis on both ergonomic working and patient comfort and accessibility; its exceptionally low minimum height makes work easier for nurses and makes getting out of bed easier for patients. The back section slide mechanism also improves patient comfort. The XS bed

also features an extremely broad Trendelenburg / anti-Trendelenburg adjustment range.

Our many years' experience in designing and manufacturing hospital beds ensures that the ScanAfia XS and XTK provides a long service life with minimal need for maintenance.

"Totally NEW model with comprehensive easy-to-use features."



ScanAfia XS Hospital Bed with side rails (short version) and electric central locking castors



ScanAfia XTK Hospital Bed with double (twin) castors, side rails (full length)

LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness

ScanAfia X-series	ScanAfia XTK-280/290	ScanAfia XTK-480/490	ScanAfia XS-280/290	ScanAfia XS-480/490
Length, cm	215	215	215	215
Width (incl. side rails)	89,5/99,5 cm	89,5/99,5 cm	89,5/99,5 cm	89,5/99,5 cm
Height adjustment (*depending on the castor, pls see matrix below)	30..82	30..82	30..82	30..82
Weight, kg	102/104	107/109	104/106	109/111
Recommended mattress size, cm	205 x 78/88 x 13	205 x 78/88 x 13	205 x 78/88 x 13	205 x 78/88 x 13
Safe working load (SWL), kg	260 kg	260 kg	260 kg	260 kg
Central locking standard Ø125mm castors	X	X	X	X
Central locking design twin Ø125mm castors	•	•	•	•
Centrally locking standard Ø150mm castors	•	•	•	•
Centrally locking design twin Ø150mm castors	•	•	•	•
Centrally locking twin Ø125mm castors, electric	-	-	•	•
Electric system, 24 V	X	X	X	X
Back section adjustment, 0° - 70°	X	X	X	X
Foot section angle, 0° - 30° (only for 4 sections models)	-	X	-	X
Knee break angle, 0° - 18° (only for 4 sections models)	-	X	-	X
Trendelenburg/ Anti-Trendelengburg, 12°/20° with gas spring	-	-	X	X
Bed lengthening of 15cm	•	•	•	•
Lifting pole	•	•	•	•
Lift support	•	•	•	•
Urine bag holder	•	•	•	•
5th castor	•	•	•	•
Bedding stand	•	•	•	•
Monitor stand	•	•	•	•
Nurse locking panel	•	•	•	•
Nurse control panel	•	•	•	•
ICU control panel	-	•	-	•
Meal tray	•	•	•	•
Handset	X	X	X	X
Manual emergency lowering of back section	X	X	X	X
Manual (quick) CPR function	-	-	•	•
Back section slide	X	X	X	X
Breathable ABS plastic base with mattress holders	X	X	X	X
Bumper rollers	X	X	X	X
Battery use	-	-	•	•
Side rails, chromed, normal or long version	•	•	•	•
Side rails, painted, normal or long version	•	•	•	•
Bed ends, in different colors	•	•	•	•
CE-maked	X	X	X	X
Warranty, 2/ 10 years	X	X	X	X
Country of manufacture	Finland	Finland	Finland	Finland

* Height adjustment range with standard 125 mm. castors	32-80 cm		32-80 cm	
Height adjustment range with standard 150 mm. castors	34-82 cm (option)		34-82 cm (option)	
Height adjustment range with design twin 125 mm. castors	30-78 cm (option)		30-78 cm (option)	
Height adjustment range with design twin 150 mm. castors	32-80 cm (option)		32-80 cm (option)	

Standard feature = X, Optional feature = •, Not available = -