

ScanAfia XS Hospital Bed

Extremely low and easy to move

Product pictured with optional accessories.



Main features:

1. Ergonomic working position thanks to impressive height adjustment range
2. High Safe Working Load (SWL) of 260 kg
3. Easy to manoeuvre thanks to its double castors (optional)

The quality ScanAfia XS hospital bed is the safe choice for hospital wards. The bed has been designed for use by healthcare professionals in hospitals, healthcare centres and similar facilities, and it fulfils the requirements of the new EN 60601-2-52:2010 standard for medical beds.

Although it can take loads up to 260 kg (SWL), the ScanAfia XS hospital bed is still light and easy to manoeuvre, thanks to the optional choice of double castors. The design places emphasis on both ergonomic working and patient comfort and accessibility; its exceptionally low minimum height makes work easier for nurses and makes getting out of bed easier for patients. The back section slide mechanism also improves patient comfort. The bed also features an extremely broad Trendelenburg / Anti-Trendelenburg adjustment range.

LOJER® For easy care


Technical data and features	
Total width (incl. siderails)	89,5 or 99,5 cm
Total length	215 cm
Height adjustment	Min. 32 cm - max. 82 cm (depending on the castor size and type)
Recommended mattress size	205 x 78/88 x 13 cm
Total weight (without rails)	80 cm wide 108 kg, 90 cm wide 110 kg
Castors	Ø125 central locking (Ø125 mm and Ø150 mm twin castors as an option)
Safe Working Load (SWL)	260 kg
Electrical system	24V
Trendelenburg/Anti-Trendelenburg	12°/ 20°
Back section adjustment	70°
Thigh section adjustment	30°
Calf section adjustment	11°
Other standard features	Hand control manual emergency lowering of back section back section slide vertical raising motion removable bed ends: white, blue, yellow or light green chrome plated or painted side rails; normal and long versions fixing for restraining belt bumper rollers fixing for lifting pole and two IV poles
CE-marked	Yes
Warranty	2 years
Country of Manufacture	Finland

Factory installed options:

- Battery for back-up
- 5th wheel to assist moving
- 125 mm twin castors
- 150 mm twin castors

Other accessories:

- Special length (extendable +15 cm) model
- Extra height side rails
- Lifting pole with triangle handle
- Lift support
- Lift support version 2
- IV pole with one hand quick adjustment
- Foldable bed end (works also as an auxiliary tray for linen)
- Nurse locking panel (allows functions to be locked to prevent patient use)
- Electric CPR
- Urine bag holder
- Reading and dining tray, on top of the side rails
- Bed end (for head end of the bed)
- Monitor tray / end
- Flexible hand controller holder



LOJER® For easy care